



GUIDE OF THE ACTIVE PRODUCTS AGAINST VIRUSES



TOGETHER, MOBILISED TO PROTECT OURSELVES AND OTHERS

MADE IN FRANCE



MEDICAL
DEVICES



SURFACES



FLOORS



SKIN
HYGIENE

Contents

Discover all the Alkapharm brand products active against enveloped viruses, such as Coronavirus, Hepatitis C, Hepatitis B, VRSh, HIV, H1N1 virus and Herpes virus (according to European standards EN 14476 and/or EN 16777 and/or EN 17111).

3


Focus on Coronavirus

	HAND HYGIENE		4
Hand disinfection by rubbing		Alkafresh Gel	4
Supports & distributors Hand hygiene			4
	SURFACES & MEDICAL DEVICES		5
Cleaning and disinfection of surfaces & medical devices		Septalkan	5
		Cidalkan	5
Cleaning & disinfection of surfaces		Surfalkan SH	5
	FLOORS		6
Cleaning & disinfection		AlkaDDH	6
		AlkaDDS	6
	MEDICAL DEVICES		6
Cleaning & pre-disinfection in dental environments		Impressiv	6
		Alkafraise	6
Cleaning & pre-disinfection		Clinalkan	7

Dear customers,

In the current situation of the Covid-19 virus pandemic, we are doing our utmost to provide you with disinfection products. Since supply is difficult and unpredictable, we will make a point of satisfying you within the limit of availability and at the fairest price.

FOCUS ON "CORONAVIRUS"



• WHAT IS CORONAVIRUS?

Coronaviruses are a large family of viruses, which cause diseases ranging from a simple cold (some seasonal viruses are coronaviruses) to more severe diseases, such as MERS

or SRAS. The virus identified to date in China is a new coronavirus. It is called **2019-nCoV**.

• WHAT HYGIENE MEASURES ARE RECOMMENDED?

How to prepare before cleaning the premises?

To limit the spread of the virus, it is important to ensure a sufficient level of disinfection in premises occupied by people or open to the public, while preserving the health of maintenance workers. Let's not forget that hand, floor and surface cleaning is the first barrier procedure to be applied routinely before disinfection or alternatively.



• Regular hand washing.

It should be noted that to be effective and limit the spread of viruses, hands must be washed following to a rigorous protocol:

- **7 steps for 15 seconds** with water and soap: palm against palm + palm against back of hand + interlaced fingers + palm against fingers + thumbs + nails + wrists.
- **Rinse your hands** for 10 seconds.
- **Dry your hands** ideally with a single-use paper towel.

• Surface cleaning.

Cleaning protocols are generally established by the hygiene manager of the establishment.

This document highlights active products on coronavirus (standards on the vaccine virus). To find all the microbiological activity of each product, go to www.sodel-sa.eu/EN and consult the technical data sheet of each product.

Use biocides carefully. Before use, read the label and product information. Before any use, make sure it is essential, especially in places frequented by the general public. Whenever possible, favour alternative methods and products with the lowest risk for human and animal health, and for the environment.
Product type 1: Disinfectants for human hygiene purposes.
Product type 2: Disinfectants and algicides not intended for direct application on humans or animals.
Product type 4: Surfaces in contact with food and feed.



HAND HYGIENE

Alkafresh gel

GLYCERINE HYDROALCOHOLIC GEL - PT1

Use on all skin types, even the most sensitive. For use on healthy skin.

Active on Norovirus, Rotavirus, enveloped viruses, such as Coronavirus, Hepatitis C, Hepatitis B, VRSh, HIV, H1N1 virus and Herpes virus: EN 14476.

- 3 mL dose (ALKE225)
- 100 mL pocket vial (ALKE118)
- 300 mL pump bottle (ALKE122)
- 500 mL pump bottle (ALKE121)
- 1 L pump bottle (ALKE120)
- 1 L airless vial (ALKE119)
- 5 L jerrycan (ALKE206)



Stainless steel support
for 1 L airless
Reference AC0026



Plastic support
for 300 and 500 mL
pump bottles
Reference AC0025



Mobile distributor
Compatible with 500 mL and 1 L
Alkafresh gel pump bottles.
Foot pedal to avoid contact.
Location for a poster intended
for the public.
Practical: gel level check window.
Lockable.
Reference ALKE233

FIND OUR ADDITIONAL PRODUCTS FOR THE CLEANING STEP ON OUR WEBSITE WWW.SODEL-SA.EU/EN



SURFACES AND MEDICAL DEVICES

Septalkan

ALCOHOL-FREE* DISINFECTANT CLEANER FOR SURFACES AND MEDICAL DEVICES - PT2/PT4

Use by wet wiping on surfaces, equipment and non-immersible invasive and non-invasive medical devices.

Active on enveloped viruses, such as Coronavirus, Hep. C, Hep. B, VRSh, HIV, H1N1, Herpes virus, etc.: EN 16777. Active on Adenovirus, Norovirus, Rotavirus, Hep. B, Hep. C and HIV: EN 14476.

- 750 mL spray (ALKE163)
- 5 L jerrycan (ALKE164)
- 100 wipes flow pack sachet - 180x200mm (ALKE165)

* Contains less than 0.1% alcohols.



CE
0459

Cidalkan

DISINFECTANT CLEANER FOR SURFACES AND MEDICAL DEVICES - PT2/PT4

Use by wet wiping on surfaces, equipment and non-immersible invasive and non-invasive medical devices.

Active on Adenovirus, Norovirus and enveloped viruses, such as Coronavirus, Hep. C, Hep. B, VRSh, HIV, H1N1, Herpes virus, etc.: EN 16777. Active on Adenovirus, Poliovirus, Norovirus, Rotavirus, Hep. B, Hep. C, HIV, H1N1 and Herpes virus: EN 14476.

Activity specific to wipes: HPV18 and HPV16 human papillomavirus in clean conditions.

- 500 mL bottle (ALKE114)
- 1 L spray (ALKE108)
- 5 L jerrycan (ALKE109)
- 120 wipes pack (130x190mm) (ALKE115)
- 200 wipes pack (216x240mm) (ALKE101)



CE
0459

Surfalkan SH

DISINFECTANT CLEANER FOR SURFACES - PT2/PT4

Use by wet wiping on surfaces and equipments.

Active on enveloped viruses, such as Coronavirus, Hep. C, Hep. B, VRSh, HIV, H1N1, Herpes virus, etc.: EN 16777. Active on Norovirus, Rotavirus, HIV, Hep. B and Hep C: EN 14476.

- 750 mL spray (ALKE203)
- 5 L jerrycan (ALKE181)



FLOORS

AlkaDDH

DISINFECTANT CLEANER FOR FLOORS & SURFACES - PT2/PT4

Use on all types of surfaces and floors by wet mopping or with a scrubber-dryer.

Active on enveloped viruses, such as Coronavirus, Hep. C, Hep. B, VRSh, HIV, H1N1 and Herpes virus: EN 14476

- 20 mL dose (ALKE207)
- 1 L dosing bottle (ALKE132)
- 5 L jerrycan (ALKE131)

AlkaDDS

SCENTED DISINFECTANT CLEANER FOR FLOORS & SURFACES - PT2

Use on all types of surfaces and floors by wet mopping or with a scrubber-dryer.

Active on enveloped viruses, such as Coronavirus, Hep. C, Hep. B, VRSh, HIV, H1N1 and Herpes virus: EN 14476.

- 1 L dosing bottle (ALKE186)



MEDICAL DEVICES

Impressiv

DENTAL IMPRESSIONS DISINFECTANT CLEANER

For dental use by spraying on impressions including Silicone, alginates, hydrocolloids and synthetic elastomers.

Active on Adenovirus, Norovirus and enveloped viruses, such as Coronavirus, Hep. B, Hep. C, VRSh, H1N1, HIV, Herpes virus, etc.: EN 16777.
Active on Rotavirus, Adenovirus, Norovirus, Hep. B, Hep. C, HIV, H1N1 and Herpes virus: EN 14476.

- 1 L spray (ALKE107)

Alkafraise

DISINFECTANT CLEANER FOR ROTATIVE INSTRUMENTS

Use by immersion, before sterilisation, on invasive and non-invasive rotating instruments in dental (drilling instruments) and podiatry environments.

Active on enveloped viruses, such as Coronavirus, Hep. B, Hep. C, VRSh, HIV, Herpes virus, etc.: EN 17111
Active on Hep. C and Herpes virus: EN 14476. Active on Hep. B: Institut Pasteur

- 2 L jerrycan (ALKE179)



MEDICAL DEVICES

Clinalkan

ALKALINE DISINFECTANT CLEANER

Use by immersion for instruments, and invasive and non-invasive medical devices, such as flexible and rigid endoscopy equipment, TEE probes, etc.

Active on Norovirus and enveloped viruses, such as Coronavirus, Hep. B, Hep. C, VRSh, H1N1, HIV, Herpes virus, etc.: EN 14476.

- 20 mL dose (ALKE200)
- 1 L dosing bottle (ALKE198)
- 5 L jerrycan (ALKE199)





MADE IN FRANCE

We are available to advise you.

Our customer service is available by phone on +33 (0)2 31 31 10 50 for any questions concerning a product or an order, or to refer you to your sales manager.

www.sodel-sa.eu/EN

Head office: Sodel - 190 Rue René Barthélemy - 14104 Lisieux Cedex - FRANCE
Tél. : + 33 (0)2 31 31 10 50 - Fax : +33 (0)2 31 31 80 60